

# Live-cell analysis of 3D spheroids: label-free & fluorescent cell health reporters


S. L. Alcantara, G. Lovell, M. Oliver, K. Patel, T. Dale, N. Holtz, E. Endsley & D. Trezise  
Essen BioScience, Welwyn Garden City, AL7 3AX UK or Ann Arbor, MI, 48108, USA

**IncuCyte**<sup>®</sup>  
by ESSEN BIOSCIENCE

## Summary & Impact

- A growing body of evidence suggests that 3D cell models yield more translational biological insight than 2D monolayers.
- Here we describe a simple kinetic live-cell imaging approach based on brightfield, in combination with fluorescent image analysis of spheroids.
- Brightfield analysis enables the monitoring of spheroid size (proliferation) and when combined with cell health markers (cytotoxicity or apoptosis) mechanisms of cell death may be elucidated.
- Expression of fluorescent proteins provides a surrogate for cell viability, where the fluorescence increases during proliferation and decreases following treatment with cytotoxic agents.
- These assays are flexible, simple and provide automated and direct measures of tumour size and health in real-time.


## Cell Health Reporters; differentiating MOA



BF or fluorescence images at 240 h post-treatment. Yellow line represents the BF mask boundary.


- SKOV3 spheroids treated with cytotoxic agent camptothecin (CMP) or the cytostatic agent cycloheximide (CHX) and loaded with IncuCyte Cytotox Green (25 nM) or IncuCyte Caspase 3/7 Green (2.5  $\mu\text{M}$ ) cell health reagents.
- Time-courses show temporal response to each agent using BF or green channels. CRCs compare the mechanism of action (cytotoxic vs cytostatic).

## Brightfield enables label free quantification of spheroids




BF images at 240 h post-treatment. Yellow line is the mask boundary.

## Brightfield analysis; a spheroid growth assay



## Continuous Live-Cell Analysis: Methodology



A flexible assay platform that sits inside a standard tissue culture incubator. IncuCyte automatically and continuously acquires and analyzes HD phase and fluorescent images of living cells cultured in microplates, dishes, or flasks.




Fast, flexible, and powerful control hub for continuous live-cell analysis comprising image acquisition, processing, and data visualization.



A suite of non-perturbing cell labeling and reporter reagents. Includes nuclear-targeted GFP and RFPs for cell counting, no-wash caspase 3/7 substrate for apoptosis, and cell kits for angiogenesis.


## 3D Spheroid assay principle




## KEY ATTRIBUTES

- Automatic image capture
- Simple mix & read
- 96- & 384-assay formats
- Fully kinetic
- Flexible assay platform
- Multiplexing


## Robust intra- and inter-plate reproducibility



- Intra-plate variability analysis; seeding density of 310 to 7,500 cells/well of 384-well plate.
- Spheroid size is seeding density-dependent.
- Profiles of A549, HT-1080 and SKOV3 cells; consistent, differential growth kinetics.
- Inter-plate variability analysis; seeding density of 2,500 cells per well.
- Reproducible BF area metrics across plates for vehicle- and CMP-treated spheroids.



## FP expression; an alternative for cell viability measurement



## FP expression; pharmacology of cell viability

